





PRESENTS:

PERCUSSIONAIRE[®] IPV[®]/CMV INTERFACING KIT A50474-2



YOUR INTRAPULMONARY PERCUSSIVE VENTILATION (IPV[®]) PERCUSSIONATOR[®] WILL SUPERIMPOSE A PERCUSSIVE DIFFUSIVE VENTILATION UPON YOUR PROGRAMMED TIDAL VOLUME DELIVERIES TO ENHANCE THE PATIENT'S LUNG RECRUITMENT.



By adding PERCUSSIVE OSCILLATION to a basic CMV convective programming with (amplitude and rate selection) any combination of a SINUSOIDAL/CMV VENTILATION can be programmed.

Your Percussionaire® Distributor can advise you how to program your convective volume pressure Ventilator to obtain an increased level of endobronchial mechanical mixing during endobronchial gas exchange.

 Percussionaire® Corporation
 1655 Glengary Bay Road, Sagle ID 83860 USA
 P.O. Box 817, Sandpoint, ID 83864 USA
 Telephone: 208.263.2549
 Fax: 208.263.0577

 Percussionaire® Corporation
 Belgium Office
 21 Broomstraat
 1820 Perk, Belgium
 Telephone: +322.751.9134
 Fax: +322.751.0923

© Percussionaire® Corp. 2009

www.percussionaire.com

US and Foreign Patents Issued and Pending.

Document ID: F-040901 DCO 10282 06-22-10