



PERCUSSIONAIRE®
CORPORATION

Advanced Cardiopulmonary Dynamics

PRESENTS:

**THE PERCUSSIONAIRE® IPV®-2C LUNG RECRUITER
MODEL F0002-C**



- ❖ **Programmable “OSCILLATORY DEMAND CONTINUOUS POSITIVE AIRWAY PRESSURE” (OD-CPAP) for reducing the work of breathing in all patient populations from neonates through pediatrics to large adults. The “clinically proven, invasive or non-invasive cardiopulmonary stabilization device”.**
- ❖ **“OSCILLATORY DEMAND CPAP” for universal Percussive Therapeutic Lung Recruitment” (TLR™)**
- ❖ **A powerful invasive or non-invasive Intrapulmonary Percussive Ventilator (IPV)® for endobronchial secretion mobilization and clearance in all patient populations.**


The IPV-2C Percussionator® is a High Frequency Percussive Ventilator (HFPV™) with:

- ❖ Selectable percussive amplitude and frequency for the mobilization and airway clearance of retained endobronchial secretions and the resolution of diffuse patchy atelectasis during mechanical airway recruitment.
- ❖ Simultaneous high volume aerosol generation for the topical delivery of water, saline bronchodilators, vasoconstrictors and certain other IV pharmaceuticals.
- ❖ An ideal non-invasive or invasive universal ventilator for acute mechanical cardio respiratory care in all patient populations from initial resuscitation to stabilization and transport.

(OD CPAP) FOR TRANSITION INTO INTRAPULMONARY PERCUSSIVE VENTILATION (IPV®)

INDEPENDENT DEMAND CPAP FOR POST VENTILATOR WEANING

Standard Percussionaire® Phasitron® Breathing Circuits can be used with the INTUITIVE IPV-2C for invasive or non-invasive ventilation.



MANUFACTURED BY:

Percussionaire® Corporation
1655 Glengary Bay Road (UPS SHIP)
P.O. Box 817 (POSTAL)
Sandpoint, Idaho 83864 U.S.A
Telephone: 208.263.2549
Fax: 208.263-0577
© Percussionaire® Corp. 2009

www.percussionaire.com

EU REPRESENTATIVE:

Percussionaire® Corporation,
Belgium Office
21 Broomstraat
1820 Perk, Belgium
Telephone: +.322.751.9134
Fax: +.322.751.0923

